MASTERMIND®

Q.F.T.

QUALITY FILTRATION TEST

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Mastermind Q.F.T. Quality Filtration Test

INTEGRATED FILTRATION TEST SYSTEM

Mastermind Q.F.T. is a laboratory instrument capable of determining both the filterability index of wine and the FOULING INDEX of water.

Purpose:

- Better comprehension of the state of wine in order to choose the best filtration system to adopt
- Result verification after wine filtration
- Check the clogging effect of utility waters
- Noticeable savings due to improved control of all filtration and filter regeneration processes

Mastermind Q.F.T. is an instrument that calculates the more significant indexes of both wine and water through fully automatic filtration weight tests under temperature control.

Wine tests:

- Filterability Index (FI)
- Modified Filterability Index (MFI)
- Maximum Volume Filterable (V MAX)

Water Test(*):

- Fouling Index (SDI)

(*) Can be performed with the W.A. Kit, sold separately

Membrane filtering processes are characterized by high efficiency in removing all contaminants; however, they often yield a limited silt retaining capacity. Therefore, the presence of solids and colloids in the product to undergo microfiltration is often the critical point of this manner of technology. In order to optimize the process, it is crucial that product filterability make its filtration worthwhile from both a technical and economic point of view.

Wine

Filterability and the rheological properties of a clear beverage can be identified through a series of indicators put forth by J. Laurenty, implemented by J.J. Descout et coll. and perfected by M. Gaillard: fouling index (FI), modified fouling index (MFI) and maximum volume filterable (Vmax). The application of these indicators to the enological world allows the direct revelation of wine's nature in order to select the best-suited filtering system. These indexes prove very useful in the appraisal of filtering system efficiency (tangential flow, cross flow, cartridge, layer, etc.).

Water

Membrane handling and use is another highly important factor. Cleaning and regeneration cycles for Juclas systems are established and validated by an in-house R&D department, which benefits from the lengthy experience gained in the management of electrolytic, tangential filtration, inverse osmosis and microfiltration (SOFOS) systems. Yet, numerous times, water used during regeneration represents the limiting factor due to bacteria, colloidal silica, iron oxide derivatives and highly clogging particles. Through an automatic test, Mastermind Q.F.T. checks the quality of water used in cellar work cycles; and defines its Silt Density Index (SDI) a.k.a FOULING INDEX, the parameter most accounted for in establishing cake layer formation. The test procedure is fully illustrated in the document drafted by the "American Standard for Testing Material" ASTM D4189 – 2007 Standard Test Method for Silt Density Index (SDI) of Water, ASTM International, West Conshohocken, PA, 2007, DOI: 10.1520/D4189-07.

MASTERMIND QFT automatically calculates these indexes of filterability, basing itself on weight detection and temperature control. The results are highly reproducible due to the extreme level of standardization that the instrument allows. Results are displayed on a PC monitor (not supplied) and stored in a database (MS Access), or they are printed out in a report (or saved in .pdf). All software involved is supplied.

The heart of the instrument is a precision scale that was specifically developed for this application, in collaboration with GIBERTINI Elettronica Srl, world leader in the production of precision weighing instrumentation. MASTERMIND QFT is the only instrument on the market that records the temperature at which the test was conducted. This feature makes it possible to perform a reliable comparison between different sets of data taken at varying moments.

Test parameters may also be modified to perform filterability index detection on particular applications (such as juice).

SYSTEM DESCRIPTION

- Pressure reduction unit for the regulation and stabilization of intake gas
- 9 kg keg with a safety valve for the test product
- Membrane gaskets with integrated temperature sensor for flow reading
- Filtered wine-collecting cylinder
- Easy coupling (Rilsan)
- Electronic high-precision scale
- USB cable for PC connection (not supplied)
- Operating software (Windows XP or Vista – or Windows 98-2000 on request)

KIT W.A. (Water Analysis) sold separately
The FOULING INDEX set includes:

- Supplementary 25 liter keg
- Specific membrane gasket, replaceable with the standard version
- Auxiliary pneumatic valve

MAIN FEATURES OF MASTERMIND QFT SOFTWARE:

- Parameter settings for tests (volumes T0-T1..., wash volume, maximum test duration)
- Scale can be calibrated
- Automatic PC port configuration

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- Test results are stored
- Recorded test temperature
- Test database with a search by field interface

FUNCTIONALITY

- Complete and simple PC interface (analysis report printouts, etc.)

The Q.F.T. (the forerunner of the current Mastermind) has been used in the finest enology laboratories in Italy since the 90s to guarantee membrane processes in accordance to the standards set forth by accreditation bodies.

All information in this form may be modified at the discretion of Juclas s.r.l.

Base theory

The filterability and the rheological characteristics of any fluid or beverage may be determined through a few parameters:

- filterability index (fouling for water) (FI);
- modified filterability index (MFI);
- maximum volume filterable (Vmax).

These values are determined, establishing the clogging of a filtering system after a certain amount of fluid (Vmax) as premise, by means of a mathematic model identified by the equation:

$$T=K \frac{V}{V_{max}-V} \quad (1)$$

Where:

T = time

V = volume of fluid filtered in time T;

Vmax = maximum volume of fluid filterable.

Formula (1) yields:

$$V_{max} = \frac{(T_2 - T_1)}{\left(\frac{T_2}{V_2}\right) - \left(\frac{T_1}{V_1}\right)} \quad (2)$$

It is standard procedure to set the values of V₁ and V₂ respectively at 200 and 400 mL, using flat 25 mm wide membranes, having an absolute porosity of 0.65 µm. The fluid examined is decarbonized wine at 20°C with filtering ΔP at a constant of 2 bar.

Volume (ml)

The first Filtering/Fouling Index (FI) is yielded from (2), (Fig. 1):

$$FI = T_2 - 2T_1 \quad (3)$$

In order to reduce operative error, a modified filtering index (MFI) was taken into account:

$$IFM = (T_3 - T_1) - 2(T_2 - T_1) \quad (4)$$

in which T₃ is the time needed to filter 600 mL of product. Operating under the above-defined conditions, FI and Vmax are put into relation by the following equation:

$$V_{max} = 400 + \frac{400 T_1}{IF} \quad (5)$$

Considering that T₁ is generally influenced by the viscosity of the device and by caking, in a negligible degree, equation (5), for a given fluid, becomes:

$$V_{max} = \frac{a}{x} + b \quad (6)$$

in which: x = FI

These considerations are confirmed for fluids characterized by low caking factors and limited solid particles or silt; such as the case of beverages destined to microfiltration. Experimental data has confirmed the hypothesized correlation (Fig.2). Nonetheless, in a few cases, demarcating situations linked to filtration problems, there are low Vmax values at low filterability indexes (FI and MFI). This occurs when T₁ is particularly influenced by caking.

Therefore, aside from the assessment of filterability indexes, it is extremely important that the exact Maximum Volume Filterable (Vmax) be determined. Such parameter makes it possible to detect specific properties of the product that are not observable

through the other indexes; thereby substantially lowering the industrial process yield (Fig.2).

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Numerous applications have confirmed that the microfiltration of a beverage is economically feasible if the following parameters are respected:

FI and MFI 10
Vmax 3000

*FERRARINI R., BUIATTI S., ZIRONI R.
Aseptic and cold packaging of beverages: proposals to increase quality and safety characteristics. Acts from the XIV Product Sector Congress, Pescara, 27-30 September 1990.

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