

Model	ESDT-3
Temperature range	30 to 40°C
Temperature accuracy	0.5°C
Temperature resolution	0.1°C
Timer range	1 min to 99 hrs, 59 min and 59 sec
Display	6 digit LED 1 digit for step no. 1 to 6 4 digit for displaying time
Sample interval	6
Bath	6 mm thick clear acrylic water bath
Bath size	31 x 16 x 26 cm
Heater	750 watt
Sweep	180° after every set time
Dimensions(L x W x H):	34.5 x 38 x 34 cm