

NEW

Julabo
THE TEMPERATURE CONTROL COMPANY

VALEGRO™

Recirculating Coolers

The new generation of chillers from JULABO
Compact – user-friendly – modern



with natural
refrigerant

ENGLISH



High level of user-friendliness and operational reliability

The ergonomically inclined OLED display and the easily accessible and illuminated level indicator ensure easy operation and enable placement under the lab bench. This saves time in the process and space in the work area. The removable venting grid warrants simple cleaning of the condenser and ensures long-term maintenance of performance and operational safety.

**adjustable
Pump**

Adjustable pump pressure

VALEGRO recirculating coolers have an infinitely adjustable pump. The flow rate can therefore be optimally adapted to the application, which enables safe working even with pressure-sensitive glass apparatus and reduces the sound level at the same time.



Optimal connectivity

RS232 and USB-C interfaces are integrated as standard, located directly under the front display and easy to access. The optional Ethernet interface enables integration into any network and full access to all functions via a PC.



	VALEGRO™ 350				VALEGRO™ 500							
Order No.	9610035				9610050							
Working temperature range °C	-20 ... +40											
Temperature stability °C	±0.3											
Cooling capacity °C (Bath fluid Ethanol) kW	+20 0.35	0 0.25	-10 0.18	-20 0.07	+20 0.5	0 0.37	-10 0.25	-20 0.1				
Pump capacity flow rate l/min	16				19							
Pump capacity pressure bar	0.6				0.85							
Digital interfaces	RS232, USB-C, Ethernet (optional)											
Menu functions	Timer-function, Infinitely adjustable pump, Record data, Ethernet can be activated retrospectively, Autostart											
Total dimensions WxDxH cm	23 x 45 x 45											

All details at:
julabo.com/products